


Comparison of Performance of Reversed Phase Packing Materials-1

Data# S-013



Column	Mightysil GP 150-4.6mmØ 5µm
Eluent	Acetonitrile/Water = 60/40
Flow Rate	1.0mL / min.
Detection	UV254nm
Temperature	40°C
Sample Vol.	5µL
Samples	1. Benzene 4.5mg/mL 2. Naphthalene 0.25mg/mL 3. Diphenyl 0.05mg/mL 4. Acenaphthene 0.01mg/mL 5. p-Terphenyl 0.25mg/mL