


Comparison of Performance of Reversed Phase Packing Materials-2

Data# S-014



Column	Mightysil GP 150-4.6mmØ 5µm
Eluent	Acetonitrile/10mM Phosphoric acid = 30/70
Flow Rate	1.0mL / min.
Detection	UV230nm
Temperature	40°C
Sample Vol.	5µL
Samples	1. Methyl-p-hydroxybenzoate 0.25mg/mL 2. Ethyl-p-hydroxybenzoate 0.25mg/mL 3. n-Propyl-p-hydroxybenzoate 0.25mg/mL 4. n-Butyl-p-hydroxybenzoate 0.25mg/mL