


Comparison of Performance of Reversed Phase Packing Materials-4

Data# S-016



Column	Mightysil GP 150-4.6mmØ 5µm
Eluent	Acetonitrile/10mM Potassium dihydrogenphosphate including 2 mM of SDS = 25/75
Flow Rate	1.0mL / min.
Detection	UV210nm
Temperature	40°C
Sample Vol.	5µL
Samples	1. Creatine phosphoric acid 0.1mg/mL 2. Creatine 0.1mg/dL 3. Creatinine 0.1mg/mL